Body Cap

Face Cap

Blackhead Cap

Body Discs

moderate (blue)

very intense (red)

Face Discs

moderate (blue)

sensitive (green)

very intense (red)

very sensitive (yellow)

Filters

GESKE

GERMAN BEAUTY TECH

Quickstart Guide


www.geske.com

GESKE Beauty Tech GmbH,
Leipziger Platz 18, 10117 Berlin, Germany
info@geske.com, www.geske.com



SmartAppGuided™
MicroDermabrasion
Blackhead Remover | 7 in 1

Short Introduction

Scan for the full instructions
in your language

Replacement Instructions

GESKE German Beauty Tech App
Available in 45 Languages

Apple and the Apple logo are trademarks of Apple Inc.
Google Play and the Google Play logo are trademarks of Google LLC.

**Read the safety information before first use.
Keep both instructions for use.**

YOUR PRODUCT AT A GLANCE

Your product is a multi-functional Microdermabrasion device. It is safe to be used up to 1 time per week, at maximum 2 times per week but not more than once per day. Please do not use more frequently. For maximal safety and ideal results, always follow your personalized Skincare Routine created by the GESKE App.

Download the GESKE app and get your personalized routine.

SPECIFICATIONS

Materials: silicone, ABS, PC, acrylic, polyester | **Storage temperature:** -5 - +60°C/23 - 140°F | **Battery:** Li-Ion 3.7V | **Battery Capacity:** 600 mAh | **Maximum Charging Current:** 370 mA | **Charging Voltage:** 5V | **Charging Time:** approx. 3 hours | **Operating Time:** Operating time differs depending on activated technologies & intensities. | **Interface:** USB | **Programs:** Spinning Exfoliation: 5 speeds, Vacuum Suction: 5 speeds, LED Light Technology: 1 colour

IMPORTANT SAFETY INFORMATION! READ BEFORE USE!

- This device has an increased scarring risk for people with Fitzpatrick Skin Types V and VI. Use at your own risk and handle with extra care.
- This device has an increased hyperpigmentation risk for people with Fitzpatrick Skin Types III-V. Use at your own risk and handle with extra care.
- Always start with the lowest intensity. A mild and steady use is desired. Higher intensity equals higher risks for unwanted effects. Only apply light pressure.
- It is highly recommended to avoid exposure to direct sunlight 3-4 weeks before any usage of the device.
- After usage of the device, direct sunlight must be avoided for the following 4-8 weeks. If you can't avoid sunlight, make sure to use sun protection with SPF 30 at minimum, SPF 50 is recommended.
- Avoid exercising and swimming and do not apply make-up for 24 h after a microdermabrasion session.
- To prevent skin infections, it is recommended to apply an antiseptic cream after use. Moisturize thoroughly.
- As a result of intensive applications with higher intensity, superficial abrasions are to be expected. The higher the intensity, the higher the risk. In case of doubts, consult your medician before any usage of the device.

The product is not suitable for use under the following conditions:

- After laser treatment or chemical peeling

- After facial surgery, both cosmetically and medically indicated

In the following cases, a consultation with a physician is necessary as to whether the product is suitable for you to use:

- If you have a heart disease
- In case of pregnancy. The safety of use during pregnancy has not been established

Further information on use:

- The product can emit intense light, avoid applying it in the direction of or near an open eye
- Never use the suction function in the same place for a longer period of time to avoid skin irritation, swelling or slight bruising
- Start with the lowest intensity and slowly increase as needed
- Do not use the product on irritated or injured skin or open wounds (after any kind of injury), in the case of bacterial or viral infections of the skin (either in the treatment area and elsewhere)
- This product is not intended to treat medical skin conditions or diseases. If you have medical skin conditions or diseases, consult your dermatologist before using this product
- Do not use the product if it is 1) untypically warm, 2) mechanically damaged or deformed, 3) discolored
- If pain and/or discomfort occurs during use, discontinue use and seek medical advice
- The eye area is very sensitive: Be extra careful in this area and avoid any contact of the product with the eyelids or the eye itself
- Only use the product as recommended in your personalized routine
- For hygiene reasons, we recommend not sharing your product with others and replacing attachments as needed
- The product is intended for cosmetic use on the face, GESKE assumes no liability for harmful consequences arising from the use on other body areas or other improper use, connection to unsuitable voltage sources, contaminated conducting liquid or balls
- Always make sure the plug and socket are completely dry before charging. Otherwise, electric shocks, a short circuit or fire may result. The product must not be used while connected to the power outlet

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